

Agicoat[®]

Agicoat Silver Calcium Alginate Dressing



Suitable for managing a wide variety of wounds including: Acute wounds, diabetic wounds and cavity and tunneling wounds with moderate to high exudate

- Accelerate the formation of granola tissue by creating optimal moisture due to the presence of carboxymethylcellulose
- Absorbs water and wound exudate up to 20 times the weight of the dressing due to the presence of Alginate calcium
- High antibacterial properties due to the presence of silver
- Soft and flexible, without sticking to the wound and so, dressing change without tissue damage.
- Due to the possibility of dryness in the wound, the dressing time on the wound is 24-48 hours.
- No need to use other disinfectants like Betadine or Antibiotic ointments

Address : Razi industrial zone, Esfahan, Iran

Phone : +98 353323396 - 9 Website: www.emadmed.com

Email: trading@emadmed.com

